

# M – FL 200 Trinocular LED fluorescence microscope

M-FL 200 is a dedicated Epi-Fluorescence microscope most suitable for beginners with affordable priced. Incorporating Infinity corrected optical system and employing plan objectives lenses as standard of image clarity and flatness in a wide range of observation methods Siedentopf observation head, double-plate stage, inward nosepiece reduces fatigue during long observations.

Upgradable with phase contrast, simple polarizing, Dark field observations with Digital Photomicrography.

# **Technical Specifications**

### **Trinocular Observation Head**

- Siedentopf observation head inclined at 30° rotatable at 360°.
- Diopter adjustment ring on ocular tube.
- Interpupillary distance from 55mm to 75mm.

#### **Eyepieces**

Extra wide field 10X (paired) F.O.V. 22mm.

## **Objectives**

# (Spherochromatic aberration infinity fluoresence)

- Flour Plan 4X
- Flour Plan 10X
- Flour Plan 40X
- Flour Plan 100X

# Mechanical Body

- Co-axial focusing system provided PRE-FOCUSING LEVER & TENSION ADJUSTMENT RING is available on right & left side respectively.
- The mechanical stage is 210mmX140mm and has co-axial low positioned controls.
- The built-in-base illumination system is incorporated.
- 3V high brightness white LED, adjustable brightness .

## **Excitation Filters**

- B Excitation (460nm-470nm) standard.
- G Excitation (520nm-530nm) standard.
- V Excitation (390nm-400nm) standard.
- UV Excitation (360nm-370nm) standard.

# **Light Source**

- Transmitted light source: 3V LED bulb.
- brightness LED lamps. Enable selecting monochromatic or dichromatic LED for light source.

#### **Standard Accessories**

Fuse tube, power cord, instruction manual, dust cover, power supply for LED lamp, blue filter and styrofoam molded pack.

# **Optional Accessories**



Digital camera attachment

Polarizer analyzer attachment

Dark field attachment

Phase contrast attachment



Medimage Technologies Pvt. Ltd. 3A, Type II, I.E., Kukatpally, Hyderabad (A.P), India - 500 072

Ph: 91-040-23075680 Fax: 91-040-23070848

Email: info@metavis.com Web: www.metavis.com

